



# AcroSeal Extra dry Solvents

Your air- and moisture-sensitive solution

# The AcroSeal system: performance, convenience, safety

## How to use AcroSeal packaging

Acros Organics™ AcroSeal™ packaging is easy to use. To get the best results, we recommend the following steps:

1. Puncture the septum
2. Withdraw the required amount
3. Replace the nitrogen blanket
4. Store upright to prevent prolonged product contact with the seal

### Advanced tips:

1. Prepare the nitrogen balloons in advance and use two balloons, one inside the other
2. Purge the doubled balloons with nitrogen three times before use
3. Bend the needle to make it easier to use and control
4. Flush each needle and syringe with nitrogen three times before use
5. Use the reaction flask to flush the syringe and needle prior to use

To get the best results from the septum, use 18- to 21-gauge needles and puncture in a clockwise pattern to avoid piercing the septum in the same place.

If you require a large amount of solvent, remove the whole cap under an inert atmosphere and decant the desired amount of product from the bottle.



## AcroSeal extra dry solvents

Choose from a wide range of extra dry solvents for moisture-sensitive and biochemistry applications. Deuterated solvents are also available in AcroSeal packaging, offering a low-cost alternative to ampule packaging and reducing the hazards associated with handling ampules.

### Available in three different grades:

- Standard grade: Suitable for most applications
- Molecular sieves grade: Stored over molecular sieves for prolonged shelf life
- Supreme grade: Filtered through 0.2 micro filters for dust-free applications

Below is the complete range of Acros Organics™ and Alfa Aesar™ extra-dry solvents in AcroSeal packaging.

Stock No.	Description	CAS No.	Brand	Package Size
32680	Acetone, 99.8%, extra dry	67-64-1	Acros Organics	1 L
32680	Acetone, 99.8%, extra dry	67-64-1	Acros Organics	100 mL
32681	Acetonitrile, 99.9%, extra dry	75-05-8	Acros Organics	1 L
32681	Acetonitrile, 99.9%, extra dry	75-05-8	Acros Organics	2.5 L
32681	Acetonitrile, 99.9%, extra dry	75-05-8	Acros Organics	100 mL
32681	Acetonitrile, 99.9%, extra dry, - FOR UCB ONLY -	75-05-8	Acros Organics	500 mL
36431	Acetonitrile, 99.9%, extra dry over molecular sieve	75-05-8	Acros Organics	1 L
36431	Acetonitrile, 99.9%, extra dry over molecular sieve	75-05-8	Acros Organics	2.5 L
36431	Acetonitrile, 99.9%, extra dry over molecular sieve	75-05-8	Acros Organics	100 mL
36431	Acetonitrile, 99.9%, extra dry over molecular sieve	75-05-8	Acros Organics	500 mL
44839	Acetonitrile, 99.9%, extra dry, package of 4 x 25 mL bottles	75-05-8	Acros Organics	100 mL
42925	Anisole, 99%, extra dry	100-66-3	Acros Organics	1 L
42925	Anisole, 99%, extra dry	100-66-3	Acros Organics	100 mL
39688	Benzyl alcohol, 98+%, extra dry	100-51-6	Acros Organics	1 L
39688	Benzyl alcohol, 98+%, extra dry	100-51-6	Acros Organics	100 mL
39721	Bis(2-methoxyethyl) ether, 99+%, extra dry	111-96-6	Acros Organics	1 L
39721	Bis(2-methoxyethyl) ether, 99+%, extra dry	111-96-6	Acros Organics	100 mL
44310	Bis(2-methoxyethyl) ether, 99+%, extra dry over molecular sieve	111-96-6	Acros Organics	1 L

Stock No.	Description	CAS No.	Brand	Package Size
44310	Bis(2-methoxyethyl) ether, 99+%, extra dry over molecular sieve	111-96-6	Acros Organics	100 mL
39896	1-Butanol, 99+%, extra dry	71-36-3	Acros Organics	1 L
39896	1-Butanol, 99+%, extra dry	71-36-3	Acros Organics	100 mL
61026	sec-Butanol, 99+%, extra dry	78-92-2	Acros Organics	1 L
61026	sec-Butanol, 99+%, extra dry	78-92-2	Acros Organics	100 mL
39695	2-Butanone, 99.5%, extra dry	78-93-3	Acros Organics	1 L
39695	2-Butanone, 99.5%, extra dry	78-93-3	Acros Organics	100 mL
42923	n-Butyl acetate, 99+%, extra dry	123-86-4	Acros Organics	1 L
42923	n-Butyl acetate, 99+%, extra dry	123-86-4	Acros Organics	100 mL
37521	tert-Butyl methyl ether, 99+%, extra dry	1634-04-4	Acros Organics	1 L
37521	tert-Butyl methyl ether, 99+%, extra dry	1634-04-4	Acros Organics	100 mL
37522	tert-Butyl methyl ether, 99+%, extra dry over molecular sieve	1634-04-4	Acros Organics	1 L
37522	tert-Butyl methyl ether, 99+%, extra dry over molecular sieve	1634-04-4	Acros Organics	100 mL
39697	Chlorobenzene, 99.8%, extra dry	108-90-7	Acros Organics	1 L
39697	Chlorobenzene, 99.8%, extra dry	108-90-7	Acros Organics	100 mL
44300	Chlorobenzene, 99.5%, extra dry over molecular sieve	108-90-7	Acros Organics	1 L
44300	Chlorobenzene, 99.5%, extra dry over molecular sieve	108-90-7	Acros Organics	100 mL
43382	1-Chlorobutane, 99.5%, anhydrous	109-69-3	Acros Organics	1 L
43382	1-Chlorobutane, 99.5%, anhydrous	109-69-3	Acros Organics	100 mL
32682	Chloroform, 99.9%, extra dry, stabilized	67-66-3	Acros Organics	1 L
32682	Chloroform, 99.9%, extra dry, stabilized	67-66-3	Acros Organics	2.5 L
32682	Chloroform, 99.9%, extra dry, stabilized	67-66-3	Acros Organics	100 mL
36432	Chloroform, 99.9%, extra dry over molecular sieve, stabilized	67-66-3	Acros Organics	1 L
36432	Chloroform, 99.9%, extra dry over molecular sieve, stabilized	67-66-3	Acros Organics	2.5 L
36432	Chloroform, 99.9%, extra dry over molecular sieve, stabilized	67-66-3	Acros Organics	100 mL
36432	Chloroform, 99.9%, extra dry over molecular sieve, stabilized	67-66-3	Acros Organics	500 mL
32683	Cyclohexane, 99.5%, extra dry	110-82-7	Acros Organics	1 L
32683	Cyclohexane, 99.5%, extra dry	110-82-7	Acros Organics	2.5 L
32683	Cyclohexane, 99.5%, extra dry	110-82-7	Acros Organics	100 mL
36466	Cyclohexane, 99.5%, extra dry over molecular sieve	110-82-7	Acros Organics	1 L
36466	Cyclohexane, 99.5%, extra dry over molecular sieve	110-82-7	Acros Organics	2.5 L
36466	Cyclohexane, 99.5%, extra dry over molecular sieve	110-82-7	Acros Organics	100 mL
36466	Cyclohexane, 99.5%, extra dry over molecular sieve	110-82-7	Acros Organics	500 mL
39768	Cyclopentane, 95+%, extra dry	287-92-3	Acros Organics	1 L
39768	Cyclopentane, 95+%, extra dry	287-92-3	Acros Organics	100 mL
39725	Cyclopentyl methyl ether, 99.5%, extra dry, stabilized	5614-37-9	Acros Organics	1 L
39725	Cyclopentyl methyl ether, 99.5%, extra dry, stabilized	5614-37-9	Acros Organics	100 mL
40617	Decahydronaphthalene, 99%, mixture of cis and trans, anhydrous	91-17-8	Acros Organics	1 L
40617	Decahydronaphthalene, 99%, mixture of cis and trans, anhydrous	91-17-8	Acros Organics	100 mL
43460	Decane, 99%, anhydrous	124-18-5	Acros Organics	1 L
43460	Decane, 99%, anhydrous	124-18-5	Acros Organics	100 mL
39691	Di-n-butyl ether, 99+%, extra dry	142-96-1	Acros Organics	1 L
39691	Di-n-butyl ether, 99+%, extra dry	142-96-1	Acros Organics	2.5 L
39691	Di-n-butyl ether, 99+%, extra dry	142-96-1	Acros Organics	100 mL
39696	1,2-Dichlorobenzene, 98+%, extra dry	95-50-1	Acros Organics	1 L
39696	1,2-Dichlorobenzene, 98+%, extra dry	95-50-1	Acros Organics	100 mL
32684	1,2-Dichloroethane, 99.8%, extra dry	107-06-2	Acros Organics	1 L
32684	1,2-Dichloroethane, 99.8%, extra dry	107-06-2	Acros Organics	2.5 L
32684	1,2-Dichloroethane, 99.8%, extra dry	107-06-2	Acros Organics	100 mL
44301	1,2-Dichloroethane, 99.5%, extra dry over molecular sieve	107-06-2	Acros Organics	1 L
44301	1,2-Dichloroethane, 99.5%, extra dry over molecular sieve	107-06-2	Acros Organics	100 mL
047352	Dichloromethane, anhydrous, 99.7+%, packaged under Argon in resealable bottles, stabilized with amylene	75-09-2	Alfa Aesar	250ml
047352	Dichloromethane, anhydrous, 99.7+%, packaged under Argon in resealable bottles, stabilized with amylene	75-09-2	Alfa Aesar	1 L
047352	Dichloromethane, anhydrous, 99.7+%, packaged under Argon in resealable bottles, stabilized with amylene	75-09-2	Alfa Aesar	4 x 1 L
32685	Dichloromethane, 99.9%, extra dry, stabilized	75-09-2	Acros Organics	1 L
32685	Dichloromethane, 99.9%, extra dry, stabilized	75-09-2	Acros Organics	2.5 L
32685	Dichloromethane, 99.9%, extra dry, stabilized	75-09-2	Acros Organics	100 mL
34846	Dichloromethane, 99.8%, extra dry over molecular sieve, stabilized	75-09-2	Acros Organics	1 L
34846	Dichloromethane, 99.8%, extra dry over molecular sieve, stabilized,	75-09-2	Acros Organics	2.5 L
34846	Dichloromethane, 99.8%, extra dry over molecular sieve, stabilized	75-09-2	Acros Organics	100 mL




Stock No.	Description	CAS No.	Brand	Package Size
34846	Dichloromethane, 99.8%, extra dry over molecular sieve, stabilized	75-09-2	Acros Organics	500 mL
44837	Dichloromethane, 99.9%, extra dry, stabilized, package of 4 × 25 mL bottles	75-09-2	Acros Organics	100 mL
43613	Diethyl carbonate, 99%, anhydrous	105-58-8	Acros Organics	1 L
43613	Diethyl carbonate, 99%, anhydrous	105-58-8	Acros Organics	100 mL
32686	Diethyl ether, 99.5%, extra dry, stabilized	60-29-7	Acros Organics	1 L
32686	Diethyl ether, 99.5%, extra dry, stabilized	60-29-7	Acros Organics	2.5 L
32686	Diethyl ether, 99.5%, extra dry, stabilized	60-29-7	Acros Organics	100 mL
36433	Diethyl ether, 99.5%, extra dry over molecular sieve, stabilized	60-29-7	Acros Organics	1 L
36433	Diethyl ether, 99.5%, extra dry over molecular sieve, stabilized	60-29-7	Acros Organics	2.5 L
36433	Diethyl ether, 99.5%, extra dry over molecular sieve, stabilized	60-29-7	Acros Organics	100 mL
36433	Diethyl ether, 99.5%, extra dry over molecular sieve, stabilized	60-29-7	Acros Organics	500 mL
44842	Diethyl ether, 99.5%, extra dry, stabilized, package of 4 × 25 mL bottles	60-29-7	Acros Organics	100 mL
37523	N,N-Dimethylacetamide, 99.5%, extra dry over molecular sieve	127-19-5	Acros Organics	1 L
37523	N,N-Dimethylacetamide, 99.5%, extra dry over molecular sieve	127-19-5	Acros Organics	100 mL
39635	N,N-Dimethylacetamide, 99.5%, extra dry	127-19-5	Acros Organics	1 L
39635	N,N-Dimethylacetamide, 99.5%, extra dry	127-19-5	Acros Organics	2.5 L
39635	N,N-Dimethylacetamide, 99.5%, extra dry	127-19-5	Acros Organics	100 mL
42873	Dimethyl carbonate, 99+%, extra dry	616-38-6	Acros Organics	1 L
42873	Dimethyl carbonate, 99+%, extra dry	616-38-6	Acros Organics	100 mL
047390	N,N-Dimethylformamide, anhydrous, 99.8%, packaged under argon	68-12-2	Alfa Aesar	4 × 100 mL
047390	N,N-Dimethylformamide, anhydrous, 99.8%, packaged under argon	68-12-2	Alfa Aesar	100 mL
047390	N,N-Dimethylformamide, anhydrous, 99.8%, packaged under argon	68-12-2	Alfa Aesar	1 L
047390	N,N-Dimethylformamide, anhydrous, 99.8%, packaged under argon	68-12-2	Alfa Aesar	4 × 1 L
32687	N,N-Dimethylformamide, 99.8%, extra dry	68-12-2	Acros Organics	1 L
32687	N,N-Dimethylformamide, 99.8%, extra dry	68-12-2	Acros Organics	2.5 L
32687	N,N-Dimethylformamide, 99.8%, extra dry	68-12-2	Acros Organics	100 mL
34843	N,N-Dimethylformamide, 99.8%, extra dry over molecular sieve	68-12-2	Acros Organics	1 L
34843	N,N-Dimethylformamide, 99.8%, extra dry over molecular sieve	68-12-2	Acros Organics	2.5 L
34843	N,N-Dimethylformamide, 99.8%, extra dry over molecular sieve	68-12-2	Acros Organics	100 mL
34843	N,N-Dimethylformamide, 99.8%, extra dry over molecular sieve	68-12-2	Acros Organics	500 mL
44838	N,N-Dimethylformamide, 99.8%, extra dry, package of 4 × 25 mL bottles	68-12-2	Acros Organics	100 mL
047330	Dimethyl sulfoxide, HPLC grade, 99.9+%, packaged under argon	67-68-5	Alfa Aesar	250 mL
047330	Dimethyl sulfoxide, HPLC grade, 99.9+%, packaged under argon	67-68-5	Alfa Aesar	1 L
047330	Dimethyl sulfoxide, HPLC grade, 99.9+%, packaged under argon	67-68-5	Alfa Aesar	4 × 1 L
047305	1,4-Dioxane, anhydrous, 99.8%, packaged under argon	123-91-1	Alfa Aesar	100 mL
047305	1,4-Dioxane, anhydrous, 99.8%, packaged under argon	123-91-1	Alfa Aesar	1 L
047305	1,4-Dioxane, anhydrous, 99.8%, packaged under argon	123-91-1	Alfa Aesar	* 4 × 1 L
32689	1,4-Dioxane, 99.8%, extra dry, stabilized	123-91-1	Acros Organics	1 L
32689	1,4-Dioxane, 99.8%, extra dry, stabilized	123-91-1	Acros Organics	2.5 L
32689	1,4-Dioxane, 99.8%, extra dry, stabilized	123-91-1	Acros Organics	100 mL
36434	1,4-Dioxane, 99.5%, extra dry over molecular sieve	123-91-1	Acros Organics	1 L
36434	1,4-Dioxane, 99.5%, extra dry over molecular sieve	123-91-1	Acros Organics	2.5 L
36434	1,4-Dioxane, 99.5%, extra dry over molecular sieve	123-91-1	Acros Organics	100 mL
36434	1,4-Dioxane, 99.5%, extra dry over molecular sieve	123-91-1	Acros Organics	500 mL
45017	1,4-Dioxane, 99.8%, extra dry, stabilized, package of 4 × 25 mL bottles	123-91-1	Acros Organics	100 mL
43156	1,3-Dioxolane, 99.8%, anhydrous, stabilized with 75 ppm BHT	646-06-0	Acros Organics	1 L
43156	1,3-Dioxolane, 99.8%, anhydrous, stabilized with 75 ppm BHT,	646-06-0	Acros Organics	100 mL
43459	n-Dodecane, 99%, anhydrous	112-40-3	Acros Organics	1 L
43459	n-Dodecane, 99%, anhydrous	112-40-3	Acros Organics	100 mL
39769	Ethanol, 99.5%, extra dry, absolute	64-17-5	Acros Organics	1 L
39769	Ethanol, 99.5%, extra dry, absolute	64-17-5	Acros Organics	2.5 L
39769	Ethanol, 99.5%, extra dry, absolute	64-17-5	Acros Organics	100 mL
32690	Ethyl acetate, 99.9%, extra dry,	141-78-6	Acros Organics	1 L
32690	Ethyl acetate, 99.9%, extra dry,	141-78-6	Acros Organics	2.5 L
32690	Ethyl acetate, 99.9%, extra dry,	141-78-6	Acros Organics	100 mL
36435	Ethyl acetate, 99.9%, extra dry over molecular sieve	141-78-6	Acros Organics	1 L
36435	Ethyl acetate, 99.9%, extra dry over molecular sieve	141-78-6	Acros Organics	2.5 L
36435	Ethyl acetate, 99.9%, extra dry over molecular sieve	141-78-6	Acros Organics	100 mL
43380	Ethylbenzene, 99.8%, anhydrous	100-41-4	Acros Organics	1 L
43380	Ethylbenzene, 99.8%, anhydrous	100-41-4	Acros Organics	100 mL
43381	Ethylene glycol, 99.8%, anhydrous	107-21-1	Acros Organics	1 L
43381	Ethylene glycol, 99.8%, anhydrous	107-21-1	Acros Organics	100 mL
44298	Ethylene glycol dimethyl ether, 99.5%, extra dry over molecular sieve	110-71-4	Acros Organics	1 L

Stock No.	Description	CAS No.	Brand	Package Size
44298	Ethylene glycol dimethyl ether, 99.5%, extra dry over molecular sieve	110-71-4	Acros Organics	100 mL
61035	Ethylene glycol dimethyl ether, 99+%, extra dry	110-71-4	Acros Organics	1 L
61035	Ethylene glycol dimethyl ether, 99+%, extra dry	110-71-4	Acros Organics	100 mL
42924	Ethyl formate, 98+%, extra dry	109-94-4	Acros Organics	100 mL
32691	n-Heptane, 99+%, extra dry	142-82-5	Acros Organics	1 L
32691	n-Heptane, 99+%, extra dry	142-82-5	Acros Organics	100 mL
36436	n-Heptane, 99+%, extra dry over molecular sieve	142-82-5	Acros Organics	1 L
36436	n-Heptane, 99+%, extra dry over molecular sieve	142-82-5	Acros Organics	2.5 L
36436	n-Heptane, 99+%, extra dry over molecular sieve	142-82-5	Acros Organics	100 mL
43614	n-Hexadecane, 99%, anhydrous	544-76-3	Acros Organics	1 L
43614	n-Hexadecane, 99%, anhydrous	544-76-3	Acros Organics	100 mL
32692	n-Hexane, 96+%, extra dry	110-54-3	Acros Organics	1 L
32692	n-Hexane, 96+%, extra dry	110-54-3	Acros Organics	2.5 L
32692	n-Hexane, 96+%, extra dry	110-54-3	Acros Organics	100 mL
36437	n-Hexane, 97%, extra dry over molecular sieve	110-54-3	Acros Organics	1 L
36437	n-Hexane, 97%, extra dry over molecular sieve	110-54-3	Acros Organics	2.5 L
36437	n-Hexane, 97%, extra dry over molecular sieve	110-54-3	Acros Organics	100 mL
36437	n-Hexane, 97%, extra dry over molecular sieve	110-54-3	Acros Organics	500 mL
61037	Hexanes, 99.9%, extra dry, mixture of isomers	92112-69-1	Acros Organics	1 L
61037	Hexanes, 99.9%, extra dry, mixture of isomers	92112-69-1	Acros Organics	100 mL
43386	Hexyl alcohol, 99%, anhydrous	111-27-3	Acros Organics	1 L
43386	Hexyl alcohol, 99%, anhydrous	111-27-3	Acros Organics	100 mL
32696	Isopropanol, 99.8%, extra dry	67-63-0	Acros Organics	1 L
32696	Isopropanol, 99.8%, extra dry	67-63-0	Acros Organics	2.5 L
32696	Isopropanol, 99.8%, extra dry	67-63-0	Acros Organics	100 mL
36440	Isopropanol, 99.5%, extra dry over molecular sieve	67-63-0	Acros Organics	1 L
36440	Isopropanol, 99.5%, extra dry over molecular sieve	67-63-0	Acros Organics	2.5 L
36440	Isopropanol, 99.5%, extra dry over molecular sieve	67-63-0	Acros Organics	100 mL
36440	Isopropanol, 99.5%, extra dry over molecular sieve	67-63-0	Acros Organics	500 mL
047339	Methanol, anhydrous, 99.9%, packaged under argon	67-56-1	Alfa Aesar	250 mL
047339	Methanol, anhydrous, 99.9%, packaged under argon	67-56-1	Alfa Aesar	1 L
32695	Methanol, 99.9%, extra dry	67-56-1	Acros Organics	1 L
32695	Methanol, 99.9%, extra dry	67-56-1	Acros Organics	2.5 L
32695	Methanol, 99.9%, extra dry	67-56-1	Acros Organics	100 mL
36439	Methanol, 99.8%, extra dry over molecular sieve	67-56-1	Acros Organics	1 L
36439	Methanol, 99.8%, extra dry over molecular sieve	67-56-1	Acros Organics	2.5 L
36439	Methanol, 99.8%, extra dry over molecular sieve	67-56-1	Acros Organics	100 mL
36439	Methanol, 99.8%, Extra dry over Molecular Sieve	67-56-1	Acros Organics	500 mL
44841	Methanol, 99.9%, extra dry, package of 4 x 25 mL bottles	67-56-1	Acros Organics	100 mL
39689	2-Methoxyethanol, 99+%, extra dry	109-86-4	Acros Organics	1 L
39689	2-Methoxyethanol, 99+%, extra dry	109-86-4	Acros Organics	100 mL
37183	Methyl acetate, 99+%, extra dry	79-20-9	Acros Organics	1 L
37183	Methyl acetate, 99+%, extra dry	79-20-9	Acros Organics	2.5 L
37183	Methyl acetate, 99+%, extra dry	79-20-9	Acros Organics	100 mL
39722	2-Methylbutane, 99+%, extra dry	78-78-4	Acros Organics	1 L
39722	2-Methylbutane, 99+%, extra dry	78-78-4	Acros Organics	100 mL
39724	Methylcyclohexane, 98+%, extra dry	108-87-2	Acros Organics	1 L
39724	Methylcyclohexane, 98+%, extra dry	108-87-2	Acros Organics	100 mL
39895	2-Methyl-1-propanol, 99+%, extra dry	78-83-1	Acros Organics	1 L
39895	2-Methyl-1-propanol, 99+%, extra dry	78-83-1	Acros Organics	100 mL
32693	1-Methyl-2-pyrrolidinone, 99.5%, extra dry	872-50-4	Acros Organics	1 L
32693	1-Methyl-2-pyrrolidinone, 99.5%, extra dry	872-50-4	Acros Organics	2.5 L
32693	1-Methyl-2-pyrrolidinone, 99.5%, extra dry	872-50-4	Acros Organics	100 mL
36438	1-Methyl-2-pyrrolidinone, 99.5%, extra dry over Molecular Sieve	872-50-4	Acros Organics	1 L
36438	1-Methyl-2-pyrrolidinone, 99.5%, extra dry over Molecular Sieve	872-50-4	Acros Organics	2.5 L
36438	1-Methyl-2-pyrrolidinone, 99.5%, extra dry over Molecular Sieve	872-50-4	Acros Organics	100 mL
36438	1-Methyl-2-pyrrolidinone, 99.5%, extra dry over Molecular Sieve	872-50-4	Acros Organics	500 mL
32688	Methyl sulfoxide, 99.7+%, extra dry	67-68-5	Acros Organics	1 L
32688	Methyl sulfoxide, 99.7+%, extra dry	67-68-5	Acros Organics	2.5 L
32688	Methyl sulfoxide, 99.7+%, extra dry	67-68-5	Acros Organics	100 mL
34844	Methyl sulfoxide, 99.7+%, extra dry over molecular sieve	67-68-5	Acros Organics	1 L
34844	Methyl sulfoxide, 99.7+%, extra dry over molecular sieve	67-68-5	Acros Organics	2.5 L
34844	Methyl sulfoxide, 99.7+%, extra dry over molecular sieve	67-68-5	Acros Organics	100 mL

Stock No.	Description	CAS No.	Brand	Package Size
34844	Methyl sulfoxide, 99.7+%, extra dry over molecular sieve,	67-68-5	Acros Organics	500 mL
39662	2-Methyltetrahydrofuran, 99+%, extra dry, stabilized	96-47-9	Acros Organics	1 L
39662	2-Methyltetrahydrofuran, 99+%, extra dry, stabilized	96-47-9	Acros Organics	2.5 L
39662	2-Methyltetrahydrofuran, 99+%, extra dry, stabilized	96-47-9	Acros Organics	100 mL
39663	2-Methyltetrahydrofuran, 99+%, extra dry over molecular sieve, stabilized	96-47-9	Acros Organics	1 L
39663	2-Methyltetrahydrofuran, 99+%, extra dry over molecular sieve, stabilized	96-47-9	Acros Organics	100 mL
39720	2-Methyltetrahydrofuran, 99+%, extra dry, stabilizer free	96-47-9	Acros Organics	1 L
39720	2-Methyltetrahydrofuran, 99+%, extra dry, stabilizer free	96-47-9	Acros Organics	100 mL
43563	n-Nonane, 99%, anhydrous	111-84-2	Acros Organics	100 mL
39690	n-Octane, 99+%, extra dry	111-65-9	Acros Organics	1 L
39690	n-Octane, 99+%, extra dry	111-65-9	Acros Organics	100 mL
43458	1-Octanol, 99%, anhydrous	111-87-5	Acros Organics	1 L
43458	1-Octanol, 99%, anhydrous	111-87-5	Acros Organics	100 mL
39723	n-Pentane, 99+%, extra dry	109-66-0	Acros Organics	1 L
39723	n-Pentane, 99+%, extra dry	109-66-0	Acros Organics	100 mL
44299	n-Pentane, 99+%, extra dry over molecular sieve	109-66-0	Acros Organics	1 L
44299	n-Pentane, 99+%, extra dry over molecular sieve	109-66-0	Acros Organics	100 mL
39692	Petroleum ether, extra dry, boiling range 40 to 60°C, water <50ppm	64742-49-0	Acros Organics	1 L
39692	Petroleum ether, extra dry, boiling range 40 to 60°C, water <50ppm	64742-49-0	Acros Organics	100 mL
39694	1-Propanol, 99.5%, extra dry	71-23-8	Acros Organics	1 L
39694	1-Propanol, 99.5%, extra dry	71-23-8	Acros Organics	100 mL
42926	Propylene carbonate, 99.5%, anhydrous	108-32-7	Acros Organics	100 mL
42926	Propylene carbonate, 99.5%, anhydrous	108-32-7	Acros Organics	800 mL
047323	Pyridine, anhydrous, 99.5+%, packaged under argon	110-86-1	Alfa Aesar	100 mL
047323	Pyridine, anhydrous, 99.5+%, packaged under argon	110-86-1	Alfa Aesar	1 L
33942	Pyridine, 99.5%, extra dry	110-86-1	Acros Organics	1 L
33942	Pyridine, 99.5%, extra dry	110-86-1	Acros Organics	2.5 L
33942	Pyridine, 99.5%, extra dry	110-86-1	Acros Organics	100 mL
36442	Pyridine, 99.5%, extra dry over molecular sieve	110-86-1	Acros Organics	1 L
36442	Pyridine, 99.5%, extra dry over molecular sieve	110-86-1	Acros Organics	2.5 L
36442	Pyridine, 99.5%, extra dry over molecular sieve	110-86-1	Acros Organics	100 mL
36442	Pyridine, 99.5%, extra dry over molecular sieve	110-86-1	Acros Organics	500 mL
32697	Tetrahydrofuran, 99.85%, extra dry, stabilized	109-99-9	Acros Organics	1 L
32697	Tetrahydrofuran, 99.85%, extra dry, stabilized	109-99-9	Acros Organics	2.5 L
32697	Tetrahydrofuran, 99.85%, extra dry, stabilized	109-99-9	Acros Organics	100 mL
34845	Tetrahydrofuran, 99.5%, extra dry over molecular sieve, stabilized	109-99-9	Acros Organics	1 L
34845	Tetrahydrofuran, 99.5%, extra dry over molecular sieve, stabilized	109-99-9	Acros Organics	2.5 L
34845	Tetrahydrofuran, 99.5%, extra dry over molecular sieve, stabilized	109-99-9	Acros Organics	100 mL
34845	Tetrahydrofuran, 99.5%, extra dry over molecular sieve, stabilized	109-99-9	Acros Organics	500 mL
44836	Tetrahydrofuran, 99.85%, extra dry, stabilized, package of 4 × 25 mL bottles	109-99-9	Acros Organics	100 mL
45707	Tetrahydrofuran, 99.9%, extra dry, unstabilized	109-99-9	Acros Organics	1 L
45707	Tetrahydrofuran, 99.9%, extra dry, unstabilized	109-99-9	Acros Organics	2.5 L
45707	Tetrahydrofuran, 99.9%, extra dry, unstabilized	109-99-9	Acros Organics	100 mL
047327	Tetrahydrofuran, anhydrous, 99.8+%, BHT-free, over molecular sieves, packed with argon	109-99-9	Alfa Aesar	100 mL
047327	Tetrahydrofuran, anhydrous, 99.8+%, BHT-free, over molecular sieves, packed with argon	109-99-9	Alfa Aesar	1 L
047382	Tetrahydrofuran, anhydrous, 99.8+%, stab. with 0.025% BHT, packaged under argon	109-99-9	Alfa Aesar	100 mL
047382	Tetrahydrofuran, anhydrous, 99.8+%, stab. with 0.025% BHT, packaged under argon	109-99-9	Alfa Aesar	1 L
047382	Tetrahydrofuran, anhydrous, 99.8+%, stab. with 0.025% BHT, packaged under argon	109-99-9	Alfa Aesar	4 × 1 L
047329	Tetrahydrofuran, spectrophotometric grade, 99.7+%, unstabilized, packaged under argon	109-99-9	Alfa Aesar	250 mL
047329	Tetrahydrofuran, spectrophotometric grade, 99.7+%, unstabilized, packaged under argon	109-99-9	Alfa Aesar	1 L
047329	Tetrahydrofuran, spectrophotometric grade, 99.7+%, unstabilized, packaged under argon	109-99-9	Alfa Aesar	*4 × 1 L
44308	1,2,3,4-Tetrahydronaphthalene, 99%, extra dry	119-64-2	Acros Organics	1 L
44308	1,2,3,4-Tetrahydronaphthalene, 99%, extra dry	119-64-2	Acros Organics	100 mL
32698	Toluene, 99.85%, extra dry	108-88-3	Acros Organics	1 L
32698	Toluene, 99.85%, extra dry	108-88-3	Acros Organics	2.5 L

Stock No.	Description	CAS No.	Brand	Package Size
32698	Toluene, 99.85%, extra dry	108-88-3	Acros Organics	100 mL
36441	Toluene, 99.85%, extra dry over molecular sieve	108-88-3	Acros Organics	1 L
36441	Toluene, 99.85%, extra dry over molecular sieve	108-88-3	Acros Organics	2.5 L
36441	Toluene, 99.85%, extra dry over molecular sieve	108-88-3	Acros Organics	100 mL
36441	Toluene, 99.85%, extra dry over molecular sieve	108-88-3	Acros Organics	500 mL
44840	Toluene, 99.85%, extra dry, package of 4 x 25 mL bottles	108-88-3	Acros Organics	100 mL
44303	Trichloroethylene, 99%, extra dry, stabilized	79-01-6	Acros Organics	1 L
44303	Trichloroethylene, 99%, extra dry, stabilized	79-01-6	Acros Organics	100 mL
42927	Triethyl orthoformate, 98%, anhydrous	122-51-0	Acros Organics	100 mL
42927	Triethyl orthoformate, 98%, anhydrous	122-51-0	Acros Organics	800 mL
42928	Trimethyl orthoformate, 99+%, anhydrous	149-73-5	Acros Organics	100 mL
42928	Trimethyl orthoformate, 99+%, anhydrous	149-73-5	Acros Organics	800 mL
32694	2,2,4-Trimethylpentane, 99.5%, extra dry	540-84-1	Acros Organics	1 L
32694	2,2,4-Trimethylpentane, 99.5%, extra dry	540-84-1	Acros Organics	100 mL
36467	2,2,4-Trimethylpentane, 99.5%, extra dry over molecular sieve	540-84-1	Acros Organics	1 L
36467	2,2,4-Trimethylpentane, 99.5%, extra dry over molecular sieve	540-84-1	Acros Organics	100 mL
44302	o-Xylene, 99%, extra dry	95-47-6	Acros Organics	1 L
44302	o-Xylene, 99%, extra dry	95-47-6	Acros Organics	100 mL
15893	m-Xylene, 99%, extra dry	108-38-3	Acros Organics	1 L
15893	m-Xylene, 99%, extra dry	108-38-3	Acros Organics	100 mL
39693	Xylenes, 98+%, extra dry, mixed isomers	1330-20-7	Acros Organics	1 L
39693	Xylenes, 98+%, extra dry, mixed isomers	1330-20-7	Acros Organics	100 mL





Browse the complete portfolio of products supporting research and manufacturing. For more information visit our websites: [acros.com](https://www.acros.com) [alfa.com](https://www.alfa.com)

For Research Use Only. Not for use in diagnostic procedures. © 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. LC\_AO\_EEM\_1020\_40

**ThermoFisher**  
SCIENTIFIC