

Data sheet

Clean-Easy PCR & Gel Purification Kit

Cat. No: AN0091(100 reactions)

Description

Clean-Easy PCR & Gel Purification Kit, developed as a 2-in-1, provides a rapid and efficient method to purify DNA and remove contaminants from reaction mixtures (e.g. PCR, digestion or labelling reactions) as well as extraction of DNA from both TAE and TBE agarose gels.

The DNA is bound to silica membrane in presence of chaotropic salts, washed and then eluted in a separate tube. The purified DNA is ready to be used in all demanding molecular biology applications, including restriction digestion, ligation, sequencing, transfection into mammalian cells and in vitro transcription.

Features

- **Simple** and **Just a few minutes** procedure.
- **Wide spectrum of size** fragments could be purified.
- DNA purified **Ready to use** in all molecular biology procedures.
- **Dual performance kit** for both PCR product clean-up and DNA purification from agarose gels.
- High DNA recovery yields.

Applications

- Fast purification of DNA from agarose gels.
- Fast purification of products from PCR amplification reactions.

Kit Components

Item	AN0091
Clean-Easy minispin columns	100
Collection tubes (2 mL)	100
Buffer PB	2X60 ml
PE Buffer*	25 ml
EB Buffer	10ml

* Ethanol (96%-100%) [not included] must be added prior to use as indicated on the label. After ethanol has been added, mark the bottle to indicate that this step has been completed.

Storage

Clean-Easy PCR & Gel Purification Kit should be stored at room temperature (15–25°C) for up to 12 months without any reduction in performance.

(Continued on reverse side)

Canvax Biotech, S.L. C/Astrónoma Cecilia Payne. Edif. Canvax. 14014 Córdoba, Spain.

☎ : +34 957 348 066

☎ : +34 957 346 217

✉ : info@canvaxbiotech.com



1- SAMPLE PREPARATION

PCR Clean Up

1. Add 5 volumes of **Buffer PB** to one volume of PCR solution and mix thoroughly by pipette.

Gel Extraction

1. Using a clean, sharp razor blade or scalpel, excise the DNA band from the agarose gel. Remove the extra agarose to reduce the size of gel slice. Place the gel slice in a 1.5 ml pre-weighted tube and weigh the gel slice (The maximum amount of gel slice per column is 400 mg).

2. Add 3 volumes of **Buffer PB** to 1 volume of gel. (For example, if the agarose gel slice is 100 mg, add 300 µl of Buffer PB)

*For gels containing more than 2% agarose, add 6 volumes of **Buffer PB** per mg of gel.*

3. Incubate at 50 °C in a water bath for 10 min or until the gel slice has completely dissolved. During incubation at 50 °C, mix by vortexing or inverting the tubes every 1 minute. **Make sure the gel slice completely dissolved. For >2% gels, increase incubation time.**

Important! For fragments <500 bp and >4 kb, add 1 volume of isopropanol to the sample and mix (For example, if the agarose gel slice is 100 mg, add 100 µl isopropanol)

2- DNA BINDING

1. Label the lid of a new spin column placed in a 2 ml collection tube. Carefully apply the mix from step 1 (1-Sample preparation) to the spin column and Centrifuge at 13000 rpm for 1 minute.
2. Place the spin column in a new 2 ml collection tube, and discard the collection tube containing the filtrate.

3- WASH

1. Add 700 µl of **buffer PE** for Wash to the minispin column and centrifuge at 13000 rpm for 1 minute.
Remember! Before using it for the first time, add ethanol (96–100%) to the PE Buffer as indicated on the bottle.
2. Discard the flow-through and centrifuge at 13000 rpm for 1 minute. *This step is essential for removing traces of PE buffer.*
3. Transfer the minispin column into a new, labelled 1.5 ml microcentrifuge tube.

4- DNA ELUTION

1. Carefully open the minispin column and pipet 30 µl **Buffer EB** or H₂O (pH=7.0-8.5) directly onto the membrane. Close the cap and incubate for 1 min at room temperature, then centrifuge at 13000 rpm for 1 min to elute DNA. *To increase the DNA yield you can warm the buffer EB/H₂O to 65 °C before adding to the column.*

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively for research purposes and *in vitro* use only. The product was not tested for use in diagnostics or for drug development, and is not suitable for administration to humans or animals. Please refer to www.canvaxbiotech.com for the Material Safety Data Sheet of the product.