

jetPRIME[®] transfection reagent

Short protocol - DNA Transfection

DAY 0: Cell seeding

→ Seed cells in **V** ml of serum containing medium according to the table below

Quantities per well, dish or flask

Culture vessel	Number of cells	V volume of serum containing medium during transfection
24-well	50 000 – 80 000	0.5 ml
6-well / 35 mm	150 000 – 250 000	2 ml
100 mm / flask 75 cm ²	1 000 000 – 2 000 000	10 ml

DAY 1: Transfection = 1:2 DNA to jetPRIME[®] reagent ratio

→ Perform transfection **in the presence of serum**

→ Use **jetPRIME[®] buffer only**


→ Transfect cells at **60-80% confluency**



Watch the video «transfection using jetPRIME[®]» on YouTube!

<http://www.youtube.com/watch?v=G39wNXaZPX4>

Dilute **X** µg of DNA
in **W** µl of jetPRIME[®] buffer
Vortex 10 s and spin down



Add transfection mix
to the cells in serum
containing medium

Vortex 10 s, spin down and
incubate 10 min at RT

If required, replace medium
4 h after transfection

Incubate 24 to 48 h

Quantities per well, dish or flask

Culture vessel	W volume of jetPRIME [®] buffer	X amount of DNA added	Y volume of jetPRIME [®] reagent
24-well	50 µl	0.5 µg	1 µl
6-well / 35 mm	200 µl	2 µg	4 µl
100 mm / flask 75 cm ²	500 µl	10 µg	20 µl

DAY 2-3: Measure gene expression

See back page for optimization tips

Download complete protocol on <http://www.polyplus-transfection.com/resources/product-literature/>

jetPRIME[®] transfection reagent

Short protocol - Optimization Tips (DNA)

+ Protocol Optimization

- + Check our online Cell Transfection Database for cell specific protocols at: <http://www.polyplus-transfection.com/resources/cell-transfection-database/>
- + Test different DNA amounts: X, 0.5X and 1.5X
- + Test different DNA/jetPRIME[®] ratios, 1:2 to 1:3



Quantities per well, dish or flask

Culture vessel	W volume of jetPRIME [®] buffer	X amount of DNA added	Y volume of jetPRIME [®] reagent
24-well	50 µl	0.25 – 0.75 µg	0.5 – 2.25 µl
6-well / 35 mm	200 µl	1 – 3 µg	2 – 9 µl
100 mm / flask 75 cm ²	500 µl	5 – 15 µg	10 – 45 µl

For HEK-293 and HeLa cells, you may decrease the DNA amount to 0.5X and use the 1:2 DNA / jetPRIME[®] ratio.

+ Tips to increase cell viability of sensitive cells

- + Replace medium after 4 h
- + Decrease DNA amount to 0.5X
- + Analyze transfection at an earlier time point (24 h after transfection instead of 48 h for instance)
- + Check that the target gene does not affect cell viability

+ Good DNA Transfection Practices

- + Store appropriately jetPRIME[®] (4°C) and the DNA
- + Ensure that cells have been passaged more than twice and less than 20 times prior to transfection. Discard overconfluent cells
- + Regularly check for mycoplasma contaminations
- + Use a reporter gene to set up and optimize transfection conditions
- + Serum quality may drastically affect transfection efficiency. When purchasing a new batch of serum or trypsin, check cell viability as well as transfection efficiency

Note: jetPRIME[®] is also recommended for virus production and DNA/siRNA cotransfection, please refer to the complete protocol available online at: <http://www.polyplus-transfection.com/resources/product-literature/>