




APPLICATION NOTE

#EP006-0



SiliaMetS[®] E-PAK[®]

Using E-PAK technology to remove chlorophyll in a biomass sample



The following study demonstrates the efficiency of E-PAK activated carbon cartridges to remove chlorophyll and other impurities from a biomass sample. Cannabis was herein chosen as the model sample.

The variables analyzed after a pre-purification step were:

- Recovery of cannabinoids (*CBD*);
- The signal level of the chlorophyll by UV-visible spectroscopy (*qualitative test*);
- The saturation of the E-PAK cartridge.

LEARN MORE

about our scavengers and E-PAK in our brochure "*Solutions for Metal & Organic Scavenging*".



Using E-PAK technology to remove chlorophyll in a biomass sample

EXPERIMENTAL METHOD


The extraction of the crude from the biomass was performed with ethanol at room temperature and at -20°C to obtain the samples for testing. Usually, extractions are done below freezing temperature to get cleaner extracts. In this experiment, an extraction at room temperature was done in order to test the efficiency of E-PAK for the removal of chlorophyll at a higher concentration.

Screening static tests (*bulk carbon suspended in ethanoic crude*) were first carried out with activated carbon at a 1:1 (*w/w*) ratio. A significant improvement in the appearance of the extract was obtained, changing from a dark cloudy solution to a clear orange color. Considering the objective was to approach large-scale conditions, subsequent tests were then carried out under dynamic conditions.

Out of the different types of carbon initially tested, SiliaCarb HA Activated Carbon gave the best results. Therefore, E-PAK SiliaCarb HA cartridges ($5 \times 1 \text{ cm}$, PN: LS-CHA-51) were chosen for the remainder of the dynamic tests. Prior to each purification process, the cartridge was conditioned with 3 column volumes, i.e. 150 mL of ethanol. Then, 350 mL of extract was passed through the cartridge. The flow rate per run was set to 7.5 mL/min to ensure optimum contact time and with a residence time of 2.5 min. Once a run was completed, the extract was collected and analyzed for its total solids, as well as the cannabinoid and chlorophyll contents (*qualitative analysis*). Using the same cartridge, this process was then repeated with new extract samples up to 8 times. The experiments were performed with solutes that were extracted at either -20°C or at room temperature to test the influence of the contaminant concentration. For extractions done on samples at -20°C , CBD was found to be less soluble getting a less concentrated solution of CBD (1.32 mg/mL) than at room temperature (2.03 mg/mL) but with lower amounts of impurities present with 0.3074 and 1.580 UV absorption, respectively.

A qualitative analysis with a UV spectrophotometer was carried out to ensure a follow-up of the chlorophyll during the passages on the E-PAK SiliaCarb HA. The maximum absorbance wavelength of chlorophyll is 665 nm in ethanol.¹ This was validated with the UV analysis of the crude extract, showing maximum absorbance at 666 nm. A picture of the different extracts was taken showing the color evolution through the different runs (*Figure 1*).

Figure 1: Successive passages of cannabis extracts (*done at room temperature*) through E-PAK cartridge.



Each cuvette corresponds to 350 ml of extract (15 g/L solutes) passed through a SiliaCarb HA Activated Carbon cartridge. From left to right: 1) extracted crude, 2) ethanol solvent, and from 3) to 10) runs 1 through 8 respectively.



RESULT

The results show that the E-PAK purification of the cannabis samples allowed to retain chlorophyll significantly. A substantial drop in the UV signal was observed compared to the crude extract.

Since the cartridge was not saturated after a single treatment (*even with highly contaminated extracts such as those obtained after extraction at room temperature*), it had been possible to easily purify several extract samples using the same E-PAK cartridge. The limiting factor for reusing cartridges is the pressure required to maintain a constant flow. It was found at the industrial scale testing that 500 L of extract could be passed through the 16.5 × 100 cm cartridge without having to increase the inlet pressure. It would also be possible to double this volume (*therefore 1000 L*), by increasing the pressure.

In respect to CBD, a very high recovery rate after each successive run was achieved, demonstrating that there is neither retention nor degradation on the cartridge. For the first extract, a recovery of only 50 % was obtained. This is explained by the fact that the activated carbon / extract ratio is very high for the first filtration, which leads to the retention of certain cannabinoids. However, the fact that the yields are higher than 100 % for the following runs demonstrates that the cannabinoids initially retained on the cartridge are released in subsequent runs. (*See Table 1*).

An average CBD yield of up to 98 % was obtained for the purification of extracts generated at room temperature although containing a higher concentration of impurities (*Table 1*). It is also consistent with the increase number of recirculations required (*greater than 6 runs*).

The end of life of the E-PAK cartridge is generally characterized by a gradual decrease and / or fluctuation in the flow. In addition, the appearance of a whitish / yellowish deposit on the carbon was observed once removed from its support (*Figure 2*).

Figure 2: E-PAK cartridge before (*left*) and after treatment (*right*).



**Table 1:** Purification results with UV (665 nm) of ethanoic crude on SiliaCarb HA (1/1) obtained at -20°C and room temperature (r. t.) with average yield after CBD purification.

| Purification results with UV (665 nm) of ethanoic crude on SiliaCarb HA (1/1) | | | | |
|---|--------------|-------|--------------------------------|--------|
| Sample | CBD (mg/mL)* | | UV Chlorophyll indice (665 nm) | |
| | -20°C | r. t. | -20°C | r. t. |
| Ethanoic crude | 1.32 | 2.03 | 0.3074 | 1.580 |
| Run 1 | 0.57 | 1.32 | 0.0022 | 0.000 |
| Run 2 | 1.60 | 2.15 | 0.0017 | 0.009 |
| Run 3 | 1.34 | 1.90 | 0.0054 | 0.0213 |
| Run 4 | 1.38 | 2.08 | 0.0092 | 0.0285 |
| Run 5 | 1.41 | 2.11 | 0.0145 | 0.0341 |
| Run 6 | 1.31 | 2.01 | 0.0183 | 0.0172 |
| Run 7 | | 1.97 | | 0.0196 |
| Run 8 | | 2.11 | | 0.0220 |
| Average yield (%) | 96 | 98 | | |

* Titration of CBD in solution. Each passage can therefore be compared to the reference value.

CONCLUSIONS

In summary, a 5 g cartridge of SiliaCarb HA Activated Carbon had the capacity of processing between 6 and 8 runs (350 mL/run) of ethanoic crude samples of CBD extracted at room temperature. When scaling up to the 16.5 × 100 cm cartridge, the treatment capacity is estimated to be 1000 L of ethanoic crude sample, which corresponds to about 200 to 250 kg of biomass and yielding up to 20 kg of purified CBD. However, this result may vary depending on the type of biomass used and the extraction conditions.

References

1. Ritchie R. J., *Photosynth. Res.* **2006**; 89; 27.

DISCOVER OUR BROCHURES

METAL & ORGANIC SCAVENGING

SiliaMetS® – Metal Scavengers
SiliaBond® – Organic Scavengers
E-PAK® – Fixed Bed Flow-Through Purification Cartridges





CHROMATOGRAPHY & PURIFICATION

SiliaFlash® – Irregular Silica Gels
SiliaSphere™ PC – Spherical Silica Gels
SiliaBond® – Chromatographic Phases
SiliaSep™ – Flash Cartridges
SiliaPlate™ – TLC Plates




SAMPLE PREPARATION

SiliaPrep™ – Silica-based SPE Cartridges & Well Plates
SiliaPrepX™ – Polymeric SPE Cartridges & Well Plates
SiliaPrep™ – Micro-SPE Tips



ANALYTICAL & PREPARATIVE CHROMATOGRAPHY

SiliaSphere™ – Spherical Silica Gels
SiliaChrom® – HPLC Columns



ORGANIC SYNTHESIS

SiliaBond® – Reagents & Oxidants



R&D SERVICES

Metal & Organic Scavenging Screening
 Organic Synthesis
 Extraction & Purification
 Custom Column Packing
 Material Science

