Plants Tissue-Clearing Reagent Suitable for Fluorescent Protein Observation iTOMEI

For observing fluorescent protein clearly without complicated procedures

Decoloring solution


Tissue-Clearing Reagent iTOMEI-D [for Plants]

5mL / 25mL [T3940]

Mounting solution

Tissue-Clearing Reagent iTOMEI-M (RI 1.40) [for Plants] 5mL / 25mL [T4003]

Improved TOMEI (iTOMEI) developed by Prof. Sakamoto *et al.* is a method that specializes TOMEI¹⁾ for fluorescent protein observation. This method enables simply transparency and clearly detecting fluorescent protein.^{2,3,4)} TCI offers suitable reagents for iTOMEI.



Advantages

- Enables transparency by simple osmosis in just a few days (more than 2 days).
- Keeps the fluorescence of a fluorescent protein such as GFP or tdTomato co-expressed as a reporter gene.
- Suppresses autofluorescence.
- Applicable to wide plant species including *Oryza sativa*, *Arabidopsis thaliana*, *Marchantia polymorpha*.






Fig. 2. Cleared *Arabidopsis thaliana* by "iTOME!" (left) Non-treatment (right) iTOMEI treatment

References

1) J. Hasegawa, Y. Sakamoto, S. Matsunaga, et al., Plant Cell Physiol. 2016, 57, 462. DOI: http://doi.org/10.1093/pcp/pcw027

2) Y. Sakamoto, S. Matsunaga, et al., Commun. Biol. **2022**, *5*, 12. DOI: https://doi.org/10.1038/s42003-021-02955-9


3) M. Sato, Y. Sakamoto, S. Matsunaga, H. Tsuji, et al., Int. J. Mol. Sci. 2022, 23, 40. DOI: https://doi.org/10.3390/ijms23010040

4) Y. Sakamoto, S. Matsunaga, Tokyo University of Science, Jpn. Kokai Tokkyo Koho 2020-026975, 2020.

Plants Tissue-Clearing Reagent Suitable for Fluorescent Protein Observation: iTOMEI

Example of use

Reagents



- 1% PFA/PBS (Preparation at time of use is recommended)*1
- Decoloring solution (Tissue-Clearing Reagent iTOMEI-D [for Plants])
- Mounting solution: Please use the suitable mounting solution for objective lens. Tissue-Clearing Reagent iTOMEI-M (RI 1.40) [for Plants] is optimized to silicone immersion objectives(RI=1.40). (lohexol solution is used in the example below.)

Fixina

Fix the sample with 1% PFA/PBS for 1 hour at room temperature.*1

(Deaerate in fixative solution using vacuum pump or syringe when the sample is above ground part.)

Washing

Remove PBS and add the decoloring solution, then let it gently shake for 24 hours with shading at room temperature.*1

Decoloring

Remove the fixation solution and add PBS, then let it rest for 5 minutes at room temperature. Repeat the same work twice.*1

Washing

Remove the decoloring solution and add PBS, then let it rest for 5 minutes at room temperature. Repeat the same work twice more.

Staining

Remove PBS and add the staining solution, then let it rest with shading at room temperature.*2

Washing

Remove the staining solution and add PBS, then let it rest for 5 minutes at room temperature. Repeat the same work twice more.

Clearing

Remove PBS and add the mounting solution, then let it gently shake for 1 hour with shading

Mounting

Mount the sample on a microscope slide with the mounting solution, seal by manicure and observe it.

asmic comp ed detailed the three

- light scatte *1 : Operation time is determined by the sample size and plant type.
 - *2 : The treatment time of DAPI staining is for 30 minutes with 5 μg/mL and the treatment time of Calcofluor White staining is for 10 minutes with 1 g/L of Calcofluor White M3R and 0.5 g/L of Evans Blue, but their adjustments are needed for the purpose.
 - *3: If you want to moderate the change in osmotic pressure, it is necessary to perform a gradual replacement with a low-concentration mounting solution.

For observing only fluorescent staining dyes

This reagent enables more rapid transparency (2-3 hours).

Tissue-Clearing Reagent TOMEI [for Plants]

100mL [T3530]

Related Products

Paraformaldehyde DAPI 2HCI [for Biochemical Research] **Acetic Acid Evans Blue Acridine Orange**

25g / 500g [P0018] 5mg [A2412] 300mL [A2035] 25g [E0197]

25g [A0132]

Ordering and **Customer Service**

TCI AMERICA

:800-423-8616 / 503-283-1681 : 888-520-1075 / 503-283-1987 E-mail: Sales-US@TCIchemicals.com

TCI EUROPE N.V.

: +32 (0)3 735 07 00 : +32 (0)3 735 07 01 E-mail: Sales-EU@TCIchemicals.com

TCI Deutschland GmbH

Tel : +49 (0)6196 64053-00 : +49 (0)6196 64053-01 E-mail: Sales-DE@TCIchemicals.com

Tokyo Chemical Industry UK Ltd.

: +44 (0)1865 78 45 60 E-mail : Sales-UK@TCIchemicals.com

TCI Chemicals (India) Pvt. Ltd.

Tel : 1800 425 7889 / 044-2262 0909 Fax : 044-2262 8902 E-mail: Sales-IN@TCIchemicals.com

梯希爱(上海)化成工业发展有限公司

:800-988-0390 / 021-67121386 : 021-6712-1385

E-mail: Sales-CN@TCIchemicals.com

TOKYO CHEMICAL INDUSTRY CO., LTD.

: +81 (0)3-5640-8878

E-mail: globalbusiness@TCIchemicals.com