

ESC Swab

ESC(EasySurface Checking) swab range includes a new generation of sampling devices for microbiological examination of surfaces, especially recommended for food and pharmaceutical industries.



Product features:

- prefilled swab device
- rayonbud
- wide range of diluents
- leak-proof closure
- transport and storage at room temperature
- conform to international standards (ISO 18593, ISO 17604)

References:

- ISO 18593:2018. Microbiology of the food chain – Horizontal method for surface sampling.
- ISO 17604:2015. Microbiology of the food chain – Carcass sampling for microbiological analysis.
- ISO 6887-1:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination. Part 1: General rules for the preparation of the initial suspension and decimal dilutions.



Ref.	Description	Packaging
85623	ESCswab BPW + Neutralizing 1 mL	100/pack
85615	ESCswab BPW + Neutralizing 4 mL	100/pack
85603	ESCswab Buffered Peptone Water 10 mL	100/pack
85604	ESCswab Buffered Peptone Water 5 mL	100/pack
85605	ESCswab D/E Neutralizing Broth 10 mL	100/pack
85606	ESCswab D/E Neutralizing Broth 5 mL	100/pack
85614	ESCswab Demi Fraser Broth 10 mL	100/pack
85621	ESCswab Lethen Broth 1 mL	100/pack
85607	ESCswab Lethen Broth 10 mL	100/pack
85608	ESCswab Lethen Broth 5 mL	100/pack
85609	ESCswab Maximum Recovery Diluent 10 mL	100/pack
85610	ESCswab Maximum Recovery Diluent 5 mL	100/pack
85612	ESCswab Mycoplasma Transport Broth 5 mL	100/pack
85622	ESCswab Neutralizing Rinse Solution 1 mL	100/pack
85601	ESCswab Neutralizing Rinse Solution 10 mL	100/pack
85602	ESCswab Neutralizing Rinse Solution 5 mL	100/pack



LIOFILCHEM's products are available in Switzerland through

CHEMIE BRUNSCHWIG AG
Auf dem Wolf 10 - 4052 Basel
sales@brunschwig-ch.com
www.chemie-brunschwig.ch

