

In today's fast-paced world, healthcare professionals are constantly seeking efficient solutions to improve patient care and streamline processes. Our self-collection range offers a compassionate and practical approach to specimen collection, making it easier for patients while supporting healthcare providers in their mission to deliver exceptional care.

When it comes to healthcare, there's no place like home. We know that, just like you, your patients want to ensure their families are well taken care of. Our self-collection solutions allow patients to feel more comfortable and empowered by providing privacy and reducing the burden of in-person visits. These options allow for greater flexibility and peace of mind for individuals and families, while also saving valuable time for healthcare professionals. It's about making healthcare simpler, more personal, and, most importantly, more focused on results that matter.

Our self-collection products offer a seamless, reliable experience for both healthcare providers and patients.

Choose the solution that's right for you:
GOOD, BETTER, OR BEST



Convenience: Collection at the patient's convenience, without compromising quality or accuracy.



Privacy: Patients can collect samples discreetly, empowering them to take charge of their health.



Time & Money-Saving: Efficient processes that reduce the need for time-consuming in-person visits.



Experience a New Level of Comfort: No swabs needed with our innovative SampleSUITE™

3 FORMATS
GOOD, BETTER & BEST



U-SWAB™



U-SWAB™
Liquid Transport



SampleSUITE™
+ U-SWAB™ Liquid Transport


U-SWAB™

SELF-COLLECTION



We ensure every sample is collected and transported in optimal condition, delivering results that advance care, discovery, and progress.

WHAT YOU NEED TO KNOW



- Consumer instructions for use available on U-SWAB.CO.UK
- Tamper Proof Seal
- Shelf-Life: 24 months

Specifications

U-SWAB™ Liquid Transport Media

Swab Length: 150mm

Swab Breakpoint: 79mm

Small Scoop Length: 79mm

Large Scoop Length: 99mm


U-SWAB™

Swab Length: 120mm

Swab Breakpoint: 80mm

Scannable & Peelable Barcode Option Available

	Range	Product Code	Description
BEST	SampleSUITE™	M200121	StoolPOD™ – Blister Pouch
	SampleSUITE™	M200122	CapillaryPOD™ – 40µl – Blister Pouch
	SampleSUITE™	M200123	SalivaPOD™ – ECO – Blister Pouch
BETTER	U-SWAB™ Liquid Transport Media	M200141	U-SWAB™ MRSA Trio – 1 Standard Plastic Foam Printed Breakpoint + 1 Blue Plastic Foam + 1 Red Plastic Foam – 12x80 Vial – 1.5ml Liquid Amies Medium
	U-SWAB™ Liquid Transport Media	M200142	U-SWAB™ MRSA Duo – 1 Standard Plastic Foam Printed Breakpoint + 1 Red Plastic Foam – 12x80 Vial – 1.5ml Liquid Amies Medium
	U-SWAB™ Liquid Transport Media	M200145	U-SWAB™ Fecal – Scoop – 12x80 Vial – 2ml Cary Blair Medium
	U-SWAB™ Liquid Transport Media	M200146	U-SWAB™ Fecal – Scoop – 16x100 Vial – 2ml Cary Blair Liquid Medium
GOOD	U-SWAB™	M200147	U-SWAB™ Fecal – Scoop – 12x80 Vial – Dry
	U-SWAB™	M200148	U-SWAB™ Fecal – Scoop – 16x100 Vial – Dry
	U-SWAB™	M200101	U-SWAB™ Vaginal – Standard Flockswab – Labelled Tube
	U-SWAB™	M200102	U-SWAB™ Vaginal – Standard Flockswab – Labelled Tube – Scannable Peel Barcode



MWE's products are available in Switzerland through

CHEMIE BRUNSCHWIG AG
Auf dem Wolf 10 - 4052 Basel
sales@brunschwig-ch.com
www.chemie-brunschwig.ch



HAVE YOU THOUGHT ABOUT:

Giving your patients the flexibility to choose the self-collection device that suits them best – from the simple U-SWAB™ Vaginal to the innovative, swab-free SampleSUITE™ option in our Best tier. Empower patients with personalised care options.