

Specialty silver salts

Many minerals found in nature contain silver as a component, (e.g. Acanthite), however most silver compounds in the laboratory are synthetically produced.

Alfa Aesar, now part of Thermo Fisher Scientific, produces an extensive range of specialized silver compounds, with over 70 products available and the capability to produce large batch sizes.

- Longstanding production experience
- High purity guaranteed
- Consistent manufacturing capabilities
- Custom sizes from 5g to 200kg
- Packaged in lightproof containers
- Bulk packing sizes are customized

The extensive range is used in a multitude of applications, ranging from medical, pharmaceutical, chemical, solar energy and many more.

Stock No.	Description	Size
11660	Silver acetate, anhydrous, 99%	25g, 100g, Bulk
11420	Silver carbonate, 99.5% (metals basis)	25g, 100g, Bulk
11421	Silver chloride, 99.9% (metals basis)	25g, 100g, 500g, Bulk
11419	Silver iodide, 99.9% (metals basis)	25g, 100g, Bulk
11414	Silver nitrate, ACS, 99.9+% (metals basis)	5g, 25g, 100g, 500g, Bulk
11407	Silver(I) oxide, 99+% (metals basis)	25g, 100g, 500g, Bulk
11415	Silver phosphate, 99%	5g, 25g, Bulk
11417	Silver sulfate, ACS, 98%	25g, 100g, Bulk
11416	Silver sulfide, 99%	5g, 25g, 100g, Bulk
11539	Silver tetrafluoroborate, 99%	10g, 50g, Bulk
88722	Silver trifluoromethanesulfonate, 98%	2g, 10g, 50g, Bulk

Full product listing is available online.

Order our products online www.alfa.com

Specialty silver salts

Principle applications of silver compounds

Chemical

Many silver salts are the chemical precursors for other silver compounds or metallic silver. Silver carbonate has been employed in the Koenigs-Knorr reaction, a substitution reaction of a glycosyl halide with an alcohol to give a glycoside. Whilst the decarboxylation of silver carboxylates leads to corresponding carboxylic acid (Hunsdiecker reaction).

Pharmaceutical

The medicinal properties of silver salts have been used in a range of pharmaceutical compositions. Among its biological uses in histology, silver nitrate will demonstrate the presence of proteins and nucleic acids.

Medical and medical devices

Silver and silver salts have been used extensively throughout history for a variety of medical purposes. Recently silver salts have been incorporated into sterile dressings.

Fuel cells and solar energy materials

Silver salts are in use in a variety of electrochemical applications, such as fuel cells, in which they exhibit ionic conductivity.

Electrodes for potentiometric measurement

For potentiometric measurement, electrochemists will use a silver/silver chloride electrode as a matter of choice.

Photographic film industry

Silver compounds have been historically used in photography. The silver products used are light sensitive and undergo fast photoinduced chemical reduction to metallic silver particles.

Glass industry

Cement and coloring agents; Stained glass can be tinted yellow with silver nitrate. Silver chloride is widely employed as a transparent cement in glass-to-metal seals, for example as a sealant in cathode ray tubes.



Order our products online www.alfa.com

ThermoFisher
SCIENTIFIC

For Research Use Only. Not for use in diagnostic procedures. © 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

PF_SpecialtySilverSalts_2016_V1 BN0405161