

Axygen® PCR Microplates, Tubes, and Sealing Options

Thermal cycler compatibility



Axygen PCR Microplates, Tubes and Sealing Options for Real-Time PCR

Cat. No.	Description	Qty/Pk	Qty/Cs
Compatible with Applied Biosystems® (7000, 7300, 7500, 7700, 7900HT, ProFlex™)			
PCR-96-AB-C	96 well clear PCR microplate	10	50
PCR-96M2-HS-C	96 well clear half skirt PCR microplate with a single notch	10	50
PCR-384M2-C	384 well PCR microplate	10	50
PCR-2CP-RT-C	0.2 mL clear 8-strip caps	125	1250
UC-500	Optically clear pressure sensitive sealing film	100	500
Compatible with Applied Biosystems® (7500 Fast, 7900 Fast, StepOne™, StepOnePlus™, ViiA™ 7 Fast, QuantStudio™ 12K)			
PCR-96-LP-AB-C	96 well clear low profile PCR microplate	25	100
PCR-2CP-RT-C	0.2 mL clear 8-strip caps for real-time PCR	125	1250
UC-500	Optically clear pressure sensitive sealing film	100	500
Compatible with BioRad®/MJ Research® (iCycler, MyiQ™, iQ5™, Chromo4™, DNA Engine Opticon® 2, CFX96™ Touch™, CFX Connect™, CFX384™ Touch™, MiniOpticon™)			
PCR-96M2-HS-C	96 well clear half skirt PCR microplate with a single notch	10	50
PCR-96-LP-FLT-C	96 well clear flat top low profile non skirt PCR microplate	25	100
PCR-96-LP-FLT-W	96 well white flat top low profile non skirt PCR microplate	25	100
PCR-96-FS-C	96 well clear full skirt PCR microplate	10	50
PCR-2CP-RT-C	0.2 mL clear 8-strip caps for real-time PCR	125	1250
UC-500	Optically clear pressure sensitive sealing film	100	500
Compatible with Eppendorf® (Mastercycler® ep realplex)			
PCR-96-FLT-C	96 well clear flat top non skirt PCR microplate	25	100
PCR-96M2-HS-C	96 well clear half skirt PCR microplate with a single notch	10	50
PCR-96-FS-C	96 well clear full skirt PCR microplate	10	50
PCR-96-FS-W	96 well white full skirt PCR microplate	10	50
PCR-2CP-RT-C	0.2 mL clear 8-strip caps for real-time PCR	125	1250
UC-500	Optically clear pressure sensitive sealing film	100	500
Compatible with Roche® (Lightcycler® 480)			
PCR-96-LC480-W	96 well white PCR microplate with sealing film	10	50
PCR-384-LC480-W	384 well white PCR microplate with sealing film	10	50
Compatible with Stratagene® (Mx4000®, Mx3000®, Mx3000P®, Mx3005P®)			
PCR-96M2-HS-C	96 well clear half skirt PCR microplate with a single notch	10	50
PCR-96-FLT-C	96 well clear flat top non skirt PCR microplate	25	100
PCR-2CP-RT-C	0.2 mL clear 8-strip caps for real-time PCR	125	1250
UC-500	Optically clear pressure sensitive sealing film	100	500
Axygen Strip Tubes (compatible with most thermal cyclers)			
PCR-0208-CP-C	0.2 mL clear thin wall PCR 8-strip tubes and clear 8-strip caps	125	1250
PCR-0208-FCP-C	0.2 mL clear thin wall PCR 8-strip tubes and clear flat 8-strip caps with writing surface	125	1250
PCR-0212-FCP-C	0.2 mL clear thin wall PCR 12-strip tubes and flat 12-strip caps	80	800
PCR-0108-LP-C	0.1 mL clear thin wall 8-strip tubes compatible with ABI Fast	125	1250
Axygen Individual Tubes (compatible with most thermal cyclers)			
PCR-02-C	0.2 mL clear flat cap thin wall PCR tubes	1000	10000
PCR-02D-C	0.2 mL clear domed cap thin wall PCR tubes	1000	10000
PCR-05-C	0.5 mL clear flat cap thin wall PCR tubes	1000	10000

Axygen® PCR Microplates, Tubes, and Sealing Options for Standard PCR

Cat. No.	Description	Qty/Pk	Qty/Cs
Compatible with Applied Biosystems® (Veriti®, 9600, 9700, 2720)			
PCR-96-AB-C	96 well clear PCR microplate	10	50
PCR-96M2-HS-C	96 well clear half skirt PCR microplate with a single notch	10	50
PCR-384M2-C	384 well PCR microplate	10	50
Compatible with Applied Biosystems (Veriti®/Fast, 9800 Fast)			
PCR-96-LP-AB-C	96 well clear low profile PCR microplate	25	100
PCR-TS	PCR tight seal thermal sealing film for all thermal cyclers	100	500
Compatible with BioRad®/MJ Research® (MyCycler™, T100™, C1000™ Touch™, S1000™, Engine Tetrad® 2, PTC-100®, PTC-200, PTC-225 Tetrad, DNA Engine Dyad®, Dyad Disciple™)			
PCR-96-FLT-C	96 well clear flat top non skirt PCR microplate	25	100
PCR-96-FS-C	96 well clear full skirt PCR microplate	10	50
PCR-96-FS-W	96 well white full skirt PCR microplate	10	50
PCR-96M2-HS-C	96 well clear PCR half skirt microplate with a single notch	10	50
PCR-384M2-C	384 well PCR microplate	10	50
PCR-TS	PCR tight seal thermal sealing film for all thermal cyclers	100	500
PCR-02-FCP-C	clear flat PCR 8-strip caps with writing surface	125	1250
Compatible with Eppendorf® (MasterCycler® ep, M384, proS, pro 384, nexus gradient, nexus, nexus eco, gradient eco, nexus GX51, nexus SX1, nexus GSX1e, nexus SX1e)			
PCR-96-FS-C	96 well clear full skirt PCR microplate	10	50
PCR-96M2-HS-C	96 well clear half skirt PCR microplate with a single notch	10	50
PCR-96-FLT-C	96 well clear flat top PCR microplate	25	100
PCR-96-LP-FLT-C	96 well clear low profile flat top PCR microplate	25	100
PCR-384M2-C	384 well PCR microplate	10	50
PCR-TS	PCR tight seal thermal sealing film for all thermal cyclers	100	500
PCR-02-FCP-C	Clear flat PCR 8-strip caps with writing surface	125	1250
Axygen Strip Tubes (compatible with most thermal cyclers)			
PCR-0208-CP-C	0.2 mL clear thin wall PCR 8-strip tubes and clear 8-strip caps	125	1250
PCR-0208-FCP-C	0.2 mL clear thin wall PCR 8-strip tubes and clear flat 8-strip caps with writing surface	125	1250
PCR-0212-FCP-C	0.2 mL clear thin wall PCR 12-strip tubes and flat 12-strip caps	80	800
PCR-0108-LP-C	0.1 mL clear thin wall 8-strip tubes compatible with ABI Fast	125	1250
Axygen Individual Tubes (compatible with most thermal cyclers)			
PCR-02-C	0.2 mL clear thin wall flat cap PCR tubes	1000	10000
PCR-02D-C	0.2 mL clear thin wall domed cap PCR tubes	1000	10000
PCR-05-C	0.5 mL clear thin wall flat cap PCR tubes	1000	10000

For additional Axigen product, technical, or distributor information, please e-mail us at axgsales@corning.com, visit our website at www.corning.com/lifesciences/axigen or call 1.800.429.9436. Outside the United States, call 510.494.8900.

CORNING | **FALCON®** **AXYGEN®** **GOSSSELIN™** **PYREX®**

Warranty/Disclaimer: Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For a listing of trademarks, visit us at www.corning.com/lifesciences/trademarks. All other trademarks in this document are the property of their respective owners.

Axygen, Inc.
A Corning Subsidiary
33210 Central Ave.
Union City, CA 94587
t 800.429.9436
t 510.494.8900
f 510.494.0700

www.corning.com/lifesciences/axigen