

Do you have the need to test multiple samples each day?

Condalab would like to make your work easier with **Ready to use culture media Bags**



Adaptable to your diluter or peristaltic pump



Save time and space



Formulations according to **ISO** and **Eur. Pharma** regulations

Trypticasein Soy Broth (TSB) EP/USP/ISO



2x5 L | CAT. 6701

3x3 L | CAT. 6704

MEDIA FILL TEST AND STERILITY TEST

Liquid culture medium used in **Sterility Test** and **Media Fill Test**, according to **Eur. Pharma** and **USP**. Used as non-selective enrichment of samples.

Buffered Peptone Water ISO



2x5 L | CAT. 6702

3x3 L | CAT. 6705

5x2 L | CAT. 6707

GENERAL DILUENT IN FOOD SAFETY TESTING

The composition of this diluent is made according to the specification of the **ISO 6579** for the detection of *Salmonella* in foods and others like **ISO Standards 6785, 6887** and **8261**.

Listeria Enrichment Broth Half Fraser ISO



2x5 L | CAT. 6703

3x3 L | CAT. 6706

5x2 L | CAT. 6708

PRIMARY ENRICHMENT FOR *LISTERIA SPP.*

Used as the primary enrichment broth according to the **ISO 11290** for the detection of *Listeria spp.* Listeria Enrichment Broth Half Fraser ISO is a modification of Fraser Broth which contains half of the concentration of nalidixic acid and acriflavine to aid in the recovery of stressed cells.