

# Corning® Gosselin™ S-Blender-1

CORNING

The Corning Gosselin S-Blender-1 is a high performance blender for microbiological sample preparation. The blender, which is suitable for a range of sample dilutions between 80 and 400 mL, ensures quick and accurate blending of any sample including food, agricultural, chemical, and cosmetic.



Digital control, keypad, and LCD display



Removable stainless steel paddles

## Easy to clean:

- ▶ Removable stainless steel paddles
- ▶ Autoclavable paddles and door
- ▶ Removable blending chamber
- ▶ Integrated waste drawer

## Adjustable speed, time, and blending strength:

- ▶ Blending capacity: 80 to 400 mL bags
- ▶ Digital control, keypad, and LCD display
- ▶ Remaining time display
- ▶ Program setting mode (9 programs)
- ▶ Adjustable timer: 10 seconds to 210 seconds or continuous
- ▶ Adjustable speed (strokes/min.): slow (450), normal (550), and fast (650)

## Other Features:

- ▶ Low noise level (51dbA)
- ▶ Protection of the motor and machinery: electronic circuit breaker, no access to the chamber during operation, door-activated start/stop
- ▶ Warranty: 3 years
- ▶ Spare parts list available
- ▶ Regional service centers available for maintenance



## Ordering Information

### Corning® Gosselin™ S-Blender-1

Cat. No.	Description	Qty/Cs
S-Blender-1	Blender	1

## Accessories

Corning has developed a full range of accessories and equipment to help provide a safe and efficient working environment in your laboratory.

### Corning Gosselin S-Diluter-1

Cat. No.	Description	Qty/Cs
S-Diluter-1	Diluter	1

### Blender Bags, without Filter

Cat. No.	Description	Dimensions (H x W) (mm)	Sterile	Packaging	Qty/Box	Qty/Cs
SM1-01	80 mL blender bag	150 x 105	SAL 10 <sup>-3</sup>	25	60 x 3	4500
SM2-02	400 mL blender bag	300 x 180	SAL 10 <sup>-3</sup>	25	20 x 5	2500
SM2-01	400 mL blender bag	300 x 180	SAL 10 <sup>-3</sup>	50	10 x 5	2500
SM3-01	3500 mL blender bag	510 x 380	SAL 10 <sup>-3</sup>	25	5 x 4	500

### Blender Bags, with Lateral or Full Filter

Cat. No.	Description	Dimensions (H x W) (mm)	Pore Size (µm)	Sterile	Packaging	Qty/Box	Qty/Cs
BBAG-03	Blender bag 400 mL, lateral filter	300 x 190	51	SAL 10 <sup>-3</sup>	25	20 x 3	1500
BBAG-04	Blender bag 400 mL, lateral filter	300 x 190	51	SAL 10 <sup>-3</sup>	50	10 x 3	1500
FBAG-03	Blender bag 400 mL, full filter	300 x 190	51	SAL 10 <sup>-3</sup>	25	20 x 3	1500
FBAG-04	Blender bag 400 mL, full filter	300 x 190	51	SAL 10 <sup>-3</sup>	50	10 x 3	1500

### Clip for Blender Bag

Cat. No.	Description	Dimensions	Packaging	Qty/Box	Qty/Cs
CLIP-01	Clip for blender bag	Total length = 230 mm; clippable part = 197 mm	50	200	1000

### Metal Closure Bags, PE

Cat. No.	Description	Dimensions (H x W) (mm)	Packaging	Qty/Cs
SBAPE7518-01	120 mL metal closure bag, with writing area	178 x 76	2 boxes of 500	1000
SBA9518-01	210 mL metal closure bag	178 x 89	2 boxes of 500	1000
SBA1123-01	450 mL metal closure bag	229 x 114	2 boxes of 500	1000
SBAPE1123-01	450 mL metal closure bag, with writing area	229 x 114	2 boxes of 500	1000
SBA1423-01	650 mL metal closure bag	229 x 140	2 boxes of 500	1000
SBAPE1523-01	650 mL metal closure bag, with writing area	229 x 140	2 boxes of 500	1000
SBA1538-01	1065 mL metal closure bag	382 x 140	2 boxes of 500	1000
SBAPE1930-01	1650 mL metal closure bag, with writing area	305 x 178	4 boxes of 250	1000

### PC Kit

Cat. No.	Description	Qty/Cs
PCKIT-02	Pipet controller, 2 PP filters with PTFE membrane (0.2 and 0.45 µm), 1 silicone pipet holder	1

For more specific information on claims, visit the Certificates page at [www.corning.com/lifesciences](http://www.corning.com/lifesciences).

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

# CORNING

**Corning Incorporated**  
Life Sciences

836 North St.  
Building 300, Suite 3401  
Tewksbury, MA 01876  
t 800.492.1110  
t 978.442.2200  
f 978.442.2476

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

#### ASIA/PACIFIC

**Australia/New Zealand**  
t 61 427286832

**China**  
t 86 21 3338 4338  
f 86 21 3338 4300

**India**  
t 91 124 4604000  
f 91 124 4604099

#### Japan

t 81 3-3586 1996  
f 81 3-3586 1291

#### Korea

t 82 2-796-9500  
f 82 2-796-9300

#### Singapore

t 65 6572-9740  
f 65 6861-2913

#### Taiwan

t 886 2-2716-0338  
f 886 2-2516-7500

#### EUROPE

CSEurope@corning.com

#### France

t 0800 916 882  
f 0800 918 636

#### Germany

t 0800 101 1153  
f 0800 101 2427

#### The Netherlands

t 020 655 79 28  
f 020 659 76 73

#### United Kingdom

t 0800 376 8660  
f 0800 279 1117

#### All Other European Countries

t +31 (0) 206 59 60 51  
f +31 (0) 206 59 76 73

#### LATIN AMERICA

grupoLA@corning.com

#### Brasil

t 55 (11) 3089-7400

#### Mexico

t (52-81) 8158-8400