



# NEW PRODUCT ANNOUNCEMENT

## NEW Flu/RSV/SARS-CoV-2 Control Panel (Inactivated Swab) NOW AVAILABLE

No one can predict what will happen when flu meets COVID-19. The similar symptomology of influenza A/B, RSV, and COVID-19 may further strain our already over-burdened healthcare system. This makes the availability of tests and controls for four major respiratory viruses, including COVID-19, crucial to the health care system to allow doctors to quickly respond with the right diagnoses and treatments.

**CAT # 8246**

### Background

The Flu/RSV/SARS-CoV-2 Control Panel (Inactivated Swab) was created to meet the coming demand of flu season + COVID-19.

### Product Details

The Flu/RSV/SARS-CoV-2 Control Panel (Inactivated Swab) is designed to monitor the extraction, amplification and detection of Influenza A, Influenza B, Respiratory Syncytial Virus, and SARS-CoV-2 molecular testing assays.

- Contains six positive controls and six negative controls
  - Positive control:
    - Influenza A (H3N2) Virus
    - Influenza A (H1N1) Virus
    - Influenza B Virus
    - Respiratory Syncytial Virus
    - SARS-CoV-2
    - A549 Cells – lung epithelial cells
  - Negative Control:
    - Coxsackievirus B1
    - A549 Cells – lung epithelial cells

### Special Considerations

- Product is IVD and CE Marked
- International restrictions may apply - please connect with your sales representative for more information

### Marketing and Promotions

Find a variety of resources for our QC panels here:

- [Promotional Literature](#)
- [Technical Publications](#)
- [Video](#)