




# High yields in suspension CHO & HEK-293

FectoPRO<sup>®</sup>


- ⊕ **Amazing protein and antibody production yields**
- ⊕ **Cost-effective Transient Gene Expression using lower plasmid DNA amount**
- ⊕ **Compatible with various mammalian expression media & cells systems**
- ⊕ **Compliant with biomanufacturing guidelines**
- ⊕ **Adapted to high density cell systems**
- ⊕ **Chemically defined and animal derived component free**

FectoPRO<sup>®</sup> gives amazing antibody yields in the ExpiCHO<sup>™</sup> system:



Transfections were performed following standard protocols. IgG production was assayed at several time after transfection. Data kindly provided by ProteoGenix SAS (France).

Superior protein yields with FectoPRO<sup>®</sup> in suspension HEK-293F cells:



Transfections were performed following standard protocols. IgG3-Fc production was measured 5 days after transfection.

FectoPRO<sup>®</sup> is compatible with various serum free media (examples below):

CHO cell system	HEK-293F cell system
FreeStyle <sup>™</sup> CHO (Gibco)	Expi293 <sup>™</sup> (Gibco)
ProCHO <sup>™</sup> 4 (Lonza <sup>®</sup> )	FreeStyle <sup>™</sup> 293 (Gibco)
FreeStyle <sup>™</sup> F17 (Gibco)	FreeStyle <sup>™</sup> F17 (Gibco)
ExpiCHO <sup>™</sup> (Gibco)	HyClone <sup>™</sup> HyCell <sup>™</sup> TransF <sup>™</sup> -H (GE Healthcare <sup>™</sup> )



Product	Ref N°	Reagent size	FectoPRO <sup>®</sup> Booster size
FectoPRO <sup>®</sup>	116-001	1 mL	1 mL
	116-010	10 mL	10 mL
	116-100	100 mL	100 mL

Contact us:  
[www.polyplus-transfection.com](http://www.polyplus-transfection.com)  
 +33 (0)3 90 40 61 80  
[support@polyplus-transfection.com](mailto:support@polyplus-transfection.com)