


Multi-Tube Jet Solvent NMR and EPR Tube Washer/Dryer



5 PTFE Coated Stainless Steel Solvent Transfer Tubes

1. Prevent tube scratches
2. Allow single hand operation
3. Fit 3 mm, 4 mm and 5 mm NMR/EPR tubes

60 mL Solvent Reservoir

Solvent Channels

Flanged Reservoir Connection

1. Retaining ring eliminates joint freeze
2. E.P.D.M. seals are resistant to most organic solvents

Calibrated Length Mark for 4", 7", 8" and 9" Tube

Clear Borosilicate Glass Body

Complete Disassembly without Tools

Easy Cleaning

Multi-Tube
Washer/Dryer

ONE
washer
for **5** mixed
tubes

fits 3 mm 4 mm 5 mm
NMR/EPR Tubes

Three Connection Styles

1. #9 Silicone Stopper (WG-1209-1)
2. 24/40 Standard Taper Joint (WG-1209-J1)
3. 29/32 ISO Taper Joint (WG-1209-J2)

Product Number	Description
WG-1209-1	Multi-tube Washer/Dryer Complete, #9 Silicone Stopper Joint
WG-1209-J1	Multi-tube Washer/Dryer Complete, 24/40 Taper Joint
WG-1209-J2	Multi-tube Washer/Dryer Complete, 29/32 Taper Joint
LG-ME1	VACUUBRAND® 75 torr Oil-Free Diaphragm Vacuum Pump



Toll Free: **1-800-220-5171**
Intl.: **856-691-3200**



Toll Free: **1-800-220-1081**
Intl.: **856-691-6206**



www.wilmad-labglass.com
cs@wilmad-labglass.com