

Designed for uniform mixing across wide range of applications

CAPP® Rondo Platform Shaker

**NEW
PRODUCT**

The CAPP Rondo Platform Shaker offers efficient mixing with a durable, maintenance-free design for extended use.

Features

- Offers multiple programmable modes, such as motion direction, speed, time, dwell time, sequence and loop operation
- Different optional attachment platforms and clamps make it possible to use almost all shapes and sizes of vessels
- The Auto Power On option ensures operation continues for the remaining time after a power failure
- 2 year warranty



**BUY 1 CAPCRP-3X & GET
CAPCRP-CS9-500 BASE PLATE FOR FREE !**

Specifications

Offer valid until 19/12/26

Speed	50-300 RPM	Platform Size	330 x 330 mm
Orbital Diameter	20 mm	Orbital Diameter	IP 21
Motion	Orbital/Linear	Programmable Features	Up to 6 programs, each with 9 sub-programs
Motor Type	Brushless DC motor	Programmable Looping	Yes
Max. Loading Capacity	7.5kg (with attachment)	Dimensions (L x W x H)	420 x 355 x 110 mm
Time Setting	1 min to 99 hours and 59 mins	Weight	10 kg
Noise Level	< 60 dB	Power Consumption	24 W
Digital Display	Yes		



CAPCRP-U



CAPCRP-D




CAPCRP-C9-250



CAPCRP-C9-500



CAPCRP-F



CAPCRP-CS9-250



CAPCRP-CS9-500



CAPCRP-CS-250

Ordering Information

Cat. No.	Description
CAPCRP-3X	CAPP Rondo Platform Shaker ø 20 mm
CAPCRP-U	Universal Platform with Adjustable Bars (Flask capacity 9 x 250 mL, 6 x 500 mL)
CAPCRP-D	Disk Attachment Base Plate with Slip-resistant Tray
CAPCRP-C9-250	Base Plate with 9 Plastic Clamps: 9 x 250 mL
CAPCRP-C9-500	Base Plate with 9 Plastic Clamps: 9 x 500 mL
CAPCRP-F	Funnel Attachment
CAPCRP-CS9-250	Base Plate with 9 Stainless Steel Clamps: 9 x 250 mL
CAPCRP-CS9-500	Base Plate with 9 Stainless Steel Clamps: 9 x 500 mL
CAPCRP-CS-100	Stainless Steel Clamp for 100 mL Flask
CAPCRP-CS-150	Stainless Steel Clamp for 150 mL Flask
CAPCRP-CS-250	Stainless Steel Clamp for 250 mL Flask
CAPCRP-CS-500	Stainless Steel Clamp for 500 mL Flask
CAPCRP-CS-1000	Stainless Steel Clamp for 1000 mL Flask

Scope of Delivery

1 pc. Platform Shaker
1 pc. Power Adapter
1 pc. Spanner
1 pc. User Manual
1 pc. Warranty Card