


# GenieColor Mycoplasma Detection Kit

Detect Mycoplasma by eye


**25% OFF**

Offer ends  
30/09/2025

## Visual Mycoplasma Detection



## Excellent Performance



**Fig.1:** GenieColor Mycoplasma Detection Kit workflow and results interpretation. If mycoplasma DNA is present, the reagent color changes from pink to yellow.

**Fig.2:** The GenieColor Mycoplasma Detection Kit was tested with different cell cultures from different sources. Results showed great compatibility with multiples cell samples.

Assay Genie and Chemie Brunschwig retain the right to amend this promotion at any time

### RAPID

1-step protocol, 60 minute assay with no special instruments

### SENSITIVE

Detect > 100 mycoplasma plasmid copies/ $\mu$ L in 1  $\mu$ L of cell culture supernatant, serum or media

### VALIDATED

Detect up to 40 mycoplasma species

### VISUAL

Identify mycoplasma contamination with a pink-to-yellow color change

CODE	PRODUCT	PACK SIZE
RGEMORV0013-25	GenieColor Mycoplasma Detection Kit	25 test
RGEMORV0013-50	GenieColor Mycoplasma Detection Kit	50 test